



## Fluoride Standard Solution, 100 mg/L, 500 mL

**Product #:** 23249  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

### Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Fluoride, 100 mg/L  $\pm 2$  mg/L  $F^-$  (NIST). 500 mL bottle

#### Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifications

Concentration:	100 mg/L $\pm 2$ mg/L $F^-$
Description:	Fluoride Standard Solution
Package Type:	Bottle
Parameter:	Fluoride
Platform:	Solution
Quantity:	500 mL
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)