



Free Chlorine Chemkey Reagents (pack of 300)

Product #:	8499300
ZAR Price:	Contact Hach
Ships within 1 week	

Free Chlorine Chemkeys, range 0.04 - 4.00 mg/L Cl_2 . For SL1000 Portable Parallel Analyser (PPA). 12 boxes of 25 (300 tests total quantity).

Specifications

Method:	10260
Parameter:	Chlorine, Free
Platform:	Chemkeys
Quantity:	pk/300
Range:	0.04 - 4.00 mg/L Cl_2