



Be Right™



Total Chlorine Chemkey Reagents (box of 25)

Product #:	9429100
ZAR Price:	Contact Hach
Ships within 1 week	

Total Chlorine Chemkeys, range 0.04 - 10 mg/L Cl_2 . For SL1000 Portable Parallel Analyser (PPA).

Specifications

Instrument:	SL1000
Method:	10260
Parameter:	Chlorine, Total
Platform:	Chemkeys
Quantity:	pk/25
Range:	0.04 - 10 mg/L Cl_2