



Beaker, polypropylene, low form, 50 mL

Product #:	108041
ZAR Price:	Contact Hach
Ships within 1 week	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifications

Body material:	Polypropylene
Capacity:	50 mL
Diameter:	51.1 mm with spout
Material:	Polypropylene, Low Form
Quantity:	1 each
Units:	1 each