



## Beaker, PP, low form, 100 mL

<b>Product #:</b>	<b>108042</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Specifications

Body material:	Polypropylene
Capacity:	100 mL
Diameter:	63.0 mm with spout
Material:	Polypropylene, Low Form
Quantity:	1 each
Units:	1 each