



## Beaker, Berzelius glass, Tall Form, without spout, 200mL

Product #:	2276145
ZAR Price:	Contact Hach
Available	

200 mL Berzelius style (tall form) beakers from Kimax, are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Specifications

Body material:	Glass
Capacity:	200 mL
Material:	Glass, Berzelius, Tall Form
Units:	1 each