



## Replacement 10 mL Buret

|                    |                |
|--------------------|----------------|
| <b>Product #:</b>  | <b>2261432</b> |
| ZAR Price:         | Contact Hach   |
| Call for ship date |                |

Glass reservoir and buret with PTFE stopcock create an easy to use workstation for standard solutions. Subdivisions of 1.0 mL and 0.1 mL.

---

### Specifications

|                   |                    |
|-------------------|--------------------|
| Body material:    | Glass              |
| Capacity:         | 10 mL              |
| Special Features: | Buret, Replacement |