



Pipette, variable volume 0.1-1.0 mL

Product #:

BBP078

ZAR Price:

Contact Hach

Ships within 1 week

Hach's variable volume single-channel, air-displacement pipettes feature robust design, comfortable ergonomics, slim-line tip ejection, and single-handed operation. Outstanding accuracy with adjustments in 5 μ L increments (0.005 mL). Ideal for use with Hach TNTplus reagents. Supplied with calibration report and maintenance tool.

Pipette tips sold separately (BBP079).

Specifications

Quantity:	1 each
What's included?:	Pipette 0.1-1.0 mL, maintenance tool, calibration certificate, basic manual.
	Tips are not included and have to be purchased separately.

What's included?

Pipette 0.1-1.0 mL, maintenance tool, calibration certificate, basic manual. Tips are not included and have to be purchased separately.