



## Pipet tips, TenSette pipet 1970001, 0.1-1.0 mL, pk/1000

<b>Product #:</b>	<b>2185628</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Perform accurate, economical pipetting with Hach's Tensette pipets. Single-channel, air-displacement pipets come in two models: 0.1-1.0 mL volumes (0.1 mL increments) or 1.0-10.0 mL (1.0 mL increments). Volume selection with a simple twist of the turret.

---

### Specifications

Quantity:	1000 /pk
-----------	----------