



Pocket Dissolved Oxygen (DO) Meter

Product #: 2968800
ZAR Price: Contact Hach
Ships within 2 weeks

Pocket DO Meter measures dissolved oxygen to determine water quality in a variety of applications including wastewater, aquariums and fish hatcheries. Simultaneously displays DO and temperature readings. Unit holds and recalls minimum and maximum readings. Features automatic temperature compensation and calibration to air. Meter is ready for use and includes a soft carrying case with belt loop, two spare probe heads, electrolyte and 4 AAA batteries. Replacement electrolyte sold separately.

Sturdy design

Heavy duty polarographic probe measures dissolved oxygen and temperature.

ATC

Automatic temperature compensation readings.

Dual display

Displays D.O. and temperature at the same time.

Water resistant

IP-67 water resistant rating.

Max/Min

Maximum and minimum readings displayed.

Specifications

Accuracy 2:	± 0.4 mg/L DO
Operating temperature range:	0 - 50 °C
Power:	4-AAA Alkaline Batteries
Quantity:	1 each
Range 2:	0 - 20.0 mg/L DO
Resolution:	1 mV
Weight:	0.176 kg