



## Test kit, total chlorine, model CN-DT, 20-2000 mg/L

**Product #:**

**2471100**

ZAR Price:  
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Kit utilizes the Hach Digital Titrator and contains 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

**Kit contains all required reagents and apparatus to perform analysis in the field**

---

### Specifications

Case Style:	F
Footnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Method:	Digital Titrator/Thiosulfate
Number of tests:	100
Parameter:	Chlorine, mid and high range - As Total Cl <sub>2</sub>
Platform:	Digital Titrator
Range:	20 - 2,000 mg/L Cl <sub>2</sub>
Smallest increment:	0.2 - 5.0
Test Kit Model:	CN-DT