



Test kit, sulfite, model SU-DT

Product #:

2063300

ZAR Price:

Contact Hach

Ships within 2 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Using the Digital Titrator, this kit contains 100 tests with a range of 10-800 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	F
Method:	Digital Titrator/Iodometric
Number of tests:	100
Parameter:	Sulfite
Platform:	Digital Titrator
Range:	4 - 800 mg/L SO_3^{2-}
Smallest Increments Steps:	0.4-4.0**
Test Kit Model:	SU-DT