



Digital titrator kit

Product #:	2270900
ZAR Price:	Contact Hach
Ships within 1 week	

Accurate, convenient titrations without the bulk, fragility, or waste of a conventional buret.

Universal Digital Titrator Kit with Digital Titrator, manual, 100-ml graduated cylinder, 125 and 250 ml Erlenmeyer flasks, demineralizer bottle, clippers, and rugged carrying case.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40 - 100 typical titrations per cartridge.

Specifications

Accessories Included:	Yes
Accuracy:	± 1 % (Uncertainty of readings is 1 digit. Most samples require more than 100 digits)
Includes:	Digital Titrator with accessories in carrying case
Weight:	0.132 kg
What's included?:	Includes Digital Titrator, manual, 10-mL graduated cylinder, 125 - and 250-mL Erlenmeyer flasks, demineralizer bottle and a carrying case.

What's included?

Includes Digital Titrator, manual, 10-mL graduated cylinder, 125 - and 250-mL Erlenmeyer flasks, demineralizer bottle and a carrying case.