



Reactive Orthophosphate (as PO₄-) Test Strips, 0-50 mg/L, 50 tests

Product #:

2757150

ZAR Price:

Contact Hach

Ships within 1 week

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Specifications

| | |
|----------------------------|--|
| Parameter: | Phosphorus, Orthophosphate (reactive) - as PO ₄ |
| Platform: | Test Strip |
| Range: | 0 - 50 ppm PO ₄ |
| Smallest Increments Steps: | Steps: 0, 5, 15, 30, 50 |
| Test Kit Model: | Test Strips |
| Units: | 50 tests |