



Incubator, Culture, Low-Profile 110 V

Product #: 2619200
ZAR Price: Contact Hach
Ships within 6-8 weeks

Hach's Low-Profile Culture Model 153 Incubator exceeds performance criteria for Standard Methods/EPA/ASTM microbiological methods of an accuracy of ± 0.1 °C @ 37 °C. Five heating elements minimise "hot spots" and contribute to a uniformity of ± 0.25 °C @ 37 °C.

Features a High Limit Thermostat and Warning Light to guard against overheating during use, stainless steel interior and shelves (2), and maintenance free heating system.

Capacity: 189 Liter

Temperature Range: Ambient 5 - 70 °C (Accuracy: ± 0.1 °C @ 20 °C)

External Dimensions (W x D x H): 762 x 762 x 813 mm (30 x 30 x 32")

High Limit Thermostat and Warning Light to guard against overheating

Easy to clean stainless steel interior and shelves

Exceeds performance criteria for Standard Methods/EPA/ASTM

Maintenance free heating system

Specifications

Power requirements (Voltage):	110 V
Quantity:	1 each
What's included?:	Includes: Incubator, 2 shelves, power cord and instructional manual.

What's included?

Includes: Incubator, 2 shelves, power cord and instructional manual.