



Colony Counter, Quebec, 220 Vac, 50 Hz

Product #:	2252102
ZAR Price:	Contact Hach
Ships within 1 week	

For counting bacterial colonies. Features 11.5 cm (4.5") lens with 1.5x magnification to eliminate parallax errors; 40-watt tungsten bulbs to illuminate dark backgrounds; Wolffhuegel guide plate and 3-wire cord. Accommodates Stewart plates. Sheet metal case. Dimensions: 25 x 28 x 27 cm (10 x 11x 10.5"). CSA certified. 220V

Illumination

Illumination is provided by 40W tungsten bulbs.

Magnification

A 113mm lens offers 1.5x magnification and is positioned to eliminate parallax errors.

Specifications

Power requirements (Voltage):	220 V
What's included?:	Includes: Colony counter, guide plate, Instructional manual and power cord.

What's included?

Includes: Colony counter, guide plate, Instructional manual and power cord.